



hearTest is a certified digital audiometer that combines the capabilities of a traditional clinical audiometer with the convenience of mobile testing. Engineered for consistent results and convenience, hearTest is your go-to solution for comprehensive hearing loss assessments.

In addition, its integration with our cloud-based data management portal ensures a seamless experience for referral, data management and reporting.

Consistent results, mobility, efficiency

Conducting hearing tests on-the-go while being assured of consistent, integrated results has never been more accessible.

We're h(ear) for you

hearX has a dedicated team of Audiological experts to assist with results interpretation*

*Please note: This feature is charged at an additional cost. Please email sales@hearxgroup.com for more info.

Use hearTest for

Clinical trials

hearTest provides reliable and precise results, essential for gathering meaningful data in clinical trials and research projects. Benefit from efficient testing that is regulatory compliant, with our tailored digital audiometry solution.

Community screening

hearTest's portability, intuitive interface and integrated data management capability makes it ideal to use in large-scale hearing screening programs.

Hearing conservation

hearTest Occ Health is an extension of hearTest, which makes compliance easy, time-efficient, and cost-effective.

Self test

hearTest's digital compatibility enables contactless hearing testing. Healthcare professionals can offer home-based testing with remote support.

Healthcare professionals

Experience an accurate, easy-to-use digital solution with adjustable protocols to ensure your patients receive personalized.

What sets hearTest apart from competitor solutions?



Clinically valid tests

Evidence-based solution, validated in peer-reviewed journals and calibrated to ISO/SANS standards.



Test reliability metrics

Smart features to ensure on-site quality control and test reliability.



Environmental noise warning

Real-time monitoring of noise to alert the user of environmental noise concerns.



Conditioning functionality

Conditioning in order to prepare patients for what to expect during their test.



Efficiency

Automated testing within minutes and pre-programmed test sequences for improved efficiency.



Extended High Frequency (EHF) testing

Determining threshold shift in the high frequency range.



Easy-to-use, adjustable protocols

Adjustable protocols and hardware options cater to diverse use cases.




Real-time, integrated data management

Patient, test and facility data are instantly consolidated on mHealth Studio Cloud, our secure online database.

And even more...

- Audiogram and degree of loss classification
- Patient signature for consent
- Downloadable PDF reports
- Narrow band masking



Hardware sets you can choose from

When you purchase your hearTest solution, you can choose between 3 hardware sets to suit your requirements.



Diagnostic hardware Set A

Samsung Tab A
Sennheiser HD 280 Pro headphones
Hardcover carry case with foam inserts



Diagnostic hardware Set B

Samsung Tab A
Sennheiser HD 280 Pro headphones v3 DAC
Hardcover carry case with foam inserts



Diagnostic hardware Set C

Samsung Tab A
RadioEar DD450 headphones v3 DAC
Hardcover carry case with foam inserts

Technical specifications

Optional settings

Test paradigm: self-test / test operator mode
Narrow band masking (tones > 40dB)

Adjustable minimum testing intensity

Set A: 10 dB HL
Set B: - 10 dB HL
Set C: - 10 dB HL

Optional settings

Tablet: IP68
Headphones: Not specified

Frequency range

Set A: 125 - 8,000 Hz
Set B: 125 - 8,000 Hz
Set C: 125 - 16,000 Hz



See hearTest in action.

Contact our team to book a free demo, and join our mission as we continue to break barriers in clinical research and hearing testing.

Shipping charges are included for local shipments in the USA. Additional charges will apply for shipments outside of the USA.

Our key clients

Te Whatu Ora
Health New Zealand
Te Toka Tumai Auckland

Children's
MINNESOTA

Lifetime
HEARING CLINIC

NATIONAL FOUNDATION FOR
Deaf & Hard of Hearing

www.hearxgroup.com | US: +1 (415) 212-5500 | ZA: +27 (0) 12 030-0268
OR email sales@hearxgroup.com

USA

2140 South DuPont Highway, Camden,
Kent County Camden, DE 19934

ZA

Ground Floor, Building 2, Ashlea Gardens
Office Park, 180 Garsfontein Road, Ashlea
Gardens, Pretoria, 0081

Please note: Shipping charges are included for local shipments in the USA. Additional charges will apply for shipments outside of the USA.

hearX® is committed to protecting your privacy and ensuring that personal information is collected and used lawfully and transparently. We use the information provided to contact you about relevant products, services and content. You may unsubscribe at any time. Standard hearX Privacy Policy applies.

FDA Registered | CE-marked | ISO 13485 Certified | SAHPRA Registered | HIPAA Compliant | POPIA Compliant | GDPR Compliant



@hearxgroup



@hearX_Group



@hearX Group



@hearxgroup

© 2016-2024 hearX IP (Pty) Ltd. Copyright in this work including its layout and the designs contained herein is reserved and the making of copies is prohibited.